

Product Datasheet

SFTA3 antibody (orb352025)

Description	SFTA3 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human Surfactant-associated protein 3 protein (1-94AA)
Target	SFTA3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P0C7M3
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.