

Product Datasheet

MOB1A antibody (Biotin) (orb352036)

Description	MOB1A antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human MOB kinase activator 1A protein (2-216AA)
Target	MOB1A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9H8S9
Expiration Date	12 months from date of receipt.