

Product Datasheet

ORF2 antibody (orb352045)

Description	ORF2 antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Hepatitis E virus genotype 1 Secreted protein ORF2 protein (20-660AA)
Target	ORF2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q68985
Expiration Date	12 months from date of receipt.