

Product Datasheet

fla antibody (orb352052)

Description	fla antibody
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Borrelia burgdorferi Flagellar filament 41 kDa core protein (1-336AA)
Target	fla
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P11089
Expiration Date	12 months from date of receipt.