

Product Datasheet

Kif1a antibody (orb352053)

Description	Kif1a antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Mouse Kinesin-like protein KIF1A protein (1-361AA)
Target	Kif1a
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P33173
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.