

Product Datasheet

PAH antibody (orb35211)

Description	Rabbit polyclonal antibody to PAH. The protein belongs to the bipterin-dependent aromatic amino aci...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	PAH
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	51862
Uniprot ID	P00439
NCBI	NP_000268.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

