



Product Datasheet

PLA2G4B antibody (orb35212)

Description	Rabbit polyclonal antibody to PLA2G4B. The protein is highly expressed in brain, heart, liver, cereb...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Synthetic Peptide
Target	PLA2G4B
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	87978
Uniprot ID	POC869
NCBI	NP_001108105.1, NP_001185517.1, NP_005081.1
Dilution Range	WB: 1:1000, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Western blot analysis of in K562 cell li...