



Product Datasheet

CSGALNACT1 antibody (orb35215)

Description	Rabbit polyclonal antibody to CSGALNACT1. The protein belongs to the chondroitin N-acetylgalactosami...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	CSGALNACT1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	61294
Uniprot ID	Q8TDX6
NCBI	NP_001123990.1 , NP_060841.5
Dilution Range	WB: 1:1000, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...