

Product Datasheet

lecA antibody (HRP) (orb352154)

Description	lecA antibody (HRP)
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Pseudomonas aeruginosa PA-I galactophilic lectin protein (2-122AA)
Target	lecA
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q05097
Expiration Date	12 months from date of receipt.