

Product Datasheet

Histone H2B antibody (orb352157)

Description	Histone H2B antibody
Species/Host	Rabbit
Reactivity	Duck
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Crocodylus niloticus Histone H2B protein (2-93AA)
Target	Histone H2B
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P14001
Expiration Date	12 months from date of receipt.