

Product Datasheet

Coat protein VP1 antibody (Biotin) (orb352164)

Description	Coat protein VP1 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human parvovirus B19 Capsid protein VP2 protein (228-781AA)
Target	Capsid protein VP1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P07299
Expiration Date	12 months from date of receipt.