

Product Datasheet

BFRF3 antibody (HRP) (orb352198)

Description	BFRF3 antibody (HRP)
Species/Host	Rabbit
Reactivity	Virus
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Epstein-Barr virus Small capsomere-interacting protein (31-176AA)
Target	SCP
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P14348
Expiration Date	12 months from date of receipt.