

Product Datasheet

Reg3b antibody (Biotin) (orb352212)

Description	Reg3b antibody (Biotin)
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Rat Regenerating islet-derived protein 3-beta protein (27-175AA)
Target	Reg3b
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P25031
Expiration Date	12 months from date of receipt.