

Product Datasheet

HSD17B14 antibody (FITC) (orb352239)

Description	HSD17B14 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human 17-beta-hydroxysteroid dehydrogenase 14 protein (1-270AA)
Target	HSD17B14
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9BPX1
Expiration Date	12 months from date of receipt.