



Product Datasheet

CASP2 antibody (orb352249)


Description	CASP2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC
Immunogen	Recombinant Human Caspase-2 protein (2-452AA)
Target	CASP2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P42575
Dilution Range	IHC-P:1:1000-2000, IF:1:200-500
Expiration Date	12 months from date of receipt.



Western blot analysis of Hela(lane 1), H...



Immunofluorescence analysis of Hela cell...



Immunohistochemical staining of human pl...