

Product Datasheet

P4HA3 antibody (HRP) (orb352278)

Description	P4HA3 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Prolyl 4-hydroxylase subunit alpha-3 protein (20-544AA)
Target	P4HA3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q7Z4N8
Expiration Date	12 months from date of receipt.