

Product Datasheet

CXorf40A antibody (orb35234)

Description	CXorf40A antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human Protein CXorf40A protein (21-158AA)
Target	CXorf40A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8TE69
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.