

Product Datasheet

Allergen Car b I antibody (Biotin) (orb352352)

Description	Allergen Car b I antibody (Biotin)
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Carpinus betulus Major pollen allergen Car b 1 isoforms 1A and 1B protein (2-160AA)
Target	Major pollen allergen Car b 1 isoforms 1A and 1B
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P38949
Expiration Date	12 months from date of receipt.