

Product Datasheet

IFNT1 antibody (Biotin) (orb352356)

Description	IFNT1 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Bovine
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Bovine Interferon tau-1 protein (24-195AA)
Target	IFNT1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P15696
Expiration Date	12 months from date of receipt.