


Product Datasheet

H2-D1 antibody (orb352369)

Description	H2-D1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Mouse H-2 class I histocompatibility antigen, D-D alpha chain protein (25-311AA)
Target	H2-D1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P01900
Expiration Date	12 months from date of receipt.



Western blot analysis of A549 whole cell...