



Product Datasheet

WDTC1 antibody (orb35237)

Description	Rabbit polyclonal antibody to WDTC1. The protein, also known asADP, DCAF9, KIAA1037, is the ortholog...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Synthetic Peptide
Target	WDTC1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	75920
Uniprot ID	Q8N5D0
NCBI	NP_055838.2 , NP_001263181.1
Dilution Range	WB: 1:1000, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Western blot analysis of in HL-60 cell l...