

Product Datasheet

H2-D1 antibody (HRP) (orb352370)

Description	H2-D1 antibody (HRP)
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Mouse H-2 class I histocompatibility antigen, D-D alpha chain protein (25-311AA)
Target	H2-D1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P01900
Expiration Date	12 months from date of receipt.