

Product Datasheet

cyp125 antibody (HRP) (orb352394)

Description	cyp125 antibody (HRP)
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Mycobacterium tuberculosis Steroid C26-monooxygenase protein (1-433AA)
Target	cyp125
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P9WPP1
Expiration Date	12 months from date of receipt.