

Product Datasheet

glaB1 antibody (FITC) (orb352399)

Description	glaB1 antibody (FITC)
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	FITC
Immunogen	Recombinant Bacteroides vulgatus Alpha-1,3-galactosidase B protein (19-592AA)
Target	glaB1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	A6KWM0
Expiration Date	12 months from date of receipt.