



Product Datasheet

POLR3A antibody (orb352433)


Description	POLR3A antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ChIP, ELISA, IHC
Immunogen	Recombinant Human DNA-directed RNA polymerase III subunit RPC1 protein (392-632AA)
Target	POLR3A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O14802
Dilution Range	IHC-P:1:20-200
Expiration Date	12 months from date of receipt.



Chromatin Immunoprecipitation
293T (1.6*...



Immunohistochemistry of paraffin-embedde...



Immunohistochemistry of paraffin-embedde...