

Product Datasheet

KIR2DS1 antibody (FITC) (orb352451)

Description	KIR2DS1 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Killer cell immunoglobulin-like receptor 2DS1 protein (22-245AA)
Target	KIR2DS1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q14954
Expiration Date	12 months from date of receipt.