

Product Datasheet

Outer capsid glycoprotein VP7 antibody (FITC) (orb352455)

Description	Outer capsid glycoprotein VP7 antibody (FITC)
Species/Host	Rabbit
Reactivity	Virus
Conjugation	FITC
Immunogen	Recombinant Rotavirus A Outer capsid glycoprotein VP7 protein (51-326AA)
Target	Outer capsid glycoprotein VP7
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P10501
Expiration Date	12 months from date of receipt.