

Product Datasheet

Reg3a antibody (Biotin) (orb352460)

Description	Reg3a antibody (Biotin)
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Rat Regenerating islet-derived protein 3-alpha protein (26-174AA)
Target	Reg3a
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P35231
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.