

Product Datasheet

Prss29 antibody (orb352461)

Description	Prss29 antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Mouse Serine protease 29 protein (18-179AA)
Target	Prss29
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q99MS4
Expiration Date	12 months from date of receipt.