


Product Datasheet

CHRNA5 antibody (orb352492)

Description	CHRNA5 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Neuronal acetylcholine receptor subunit alpha-5 protein (23-240AA)
Target	CHRNA5
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P30532
Dilution Range	WB:1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat(Lane 1),...