



Product Datasheet

PIP4K2B antibody (orb352523)

Description	PIP4K2B antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human Phosphatidylinositol 5-phosphate 4-kinase type-2 beta protein (1-270AA)
Target	PIP4K2B
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P78356
Dilution Range	WB:1:500-2000, IHC-P:1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa(Lane 1), M...



Immunohistochemical staining of human br...