



Product Datasheet

TERF2IP antibody (orb352587)


Description	TERF2IP antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Recombinant Human Telomeric repeat-binding factor 2-interacting protein 1 protein (1-280AA)
Target	TERF2IP
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9NYB0
Dilution Range	IHC-P:1:20-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human bladder tissue.



Immunohistochemical staining of human brain tissue.



Immunohistochemical staining of human brain tissue.