



Product Datasheet

THOC1 antibody (orb352616)


Description	THOC1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ChIP, ELISA, IHC, WB
Immunogen	Recombinant Human THO complex subunit 1 protein (428-657AA)
Target	THOC1
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96FV9
Dilution Range	WB:1:500-2000, IHC-P:1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat cells (Lane 1)...



Immunohistochemical staining of human tissue...



Immunoprecipitation analysis of 293T cells...