


## Product Datasheet

### TNNC2 antibody (orb352621)

|                            |                                                                 |
|----------------------------|-----------------------------------------------------------------|
| <b>Description</b>         | TNNC2 antibody                                                  |
| <b>Species/Host</b>        | Rabbit                                                          |
| <b>Reactivity</b>          | Human, Mouse                                                    |
| <b>Conjugation</b>         | Unconjugated                                                    |
| <b>Tested Applications</b> | ELISA, WB                                                       |
| <b>Immunogen</b>           | Recombinant Human Troponin C, skeletal muscle protein (1-160AA) |
| <b>Target</b>              | TNNC2                                                           |
| <b>Preservatives</b>       | PBS with 0.02% sodium azide, 50% glycerol, pH7.3.               |
| <b>Form/Appearance</b>     | Liquid                                                          |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.   |
| <b>Note</b>                | For research use only                                           |
| <b>Isotype</b>             | IgG                                                             |
| <b>Clonality</b>           | Polyclonal                                                      |
| <b>Uniprot ID</b>          | <b>P02585</b>                                                   |
| <b>Dilution Range</b>      | WB:1:500-2000                                                   |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                 |



Western blot analysis of mouse skeletal ...