


Product Datasheet

TOMM70A antibody (orb352624)

Description	TOMM70A antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Mitochondrial import receptor subunit TOM70 protein (60-340AA)
Target	TOMM70
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O94826
Dilution Range	WB:1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of MCF-7 whole cel...