

Product Datasheet

UBE2D4 antibody (orb352629)

Description	UBE2D4 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human Ubiquitin-conjugating enzyme E2 D4 protein (1-147AA)
Target	UBE2D4
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9Y2X8
Expiration Date	12 months from date of receipt.