


Product Datasheet

TSEN2 antibody (orb352649)

Description	TSEN2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human tRNA-splicing endonuclease subunit Sen2 protein (1-220AA)
Target	TSEN2
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8NCE0
Dilution Range	WB:1:200-1000
Expiration Date	12 months from date of receipt.



Western blot analysis of 293T whole cell...