



Product Datasheet

ATAD3B antibody (orb352651)

Description	ATAD3B antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human ATPase family AAA domain-containing protein 3B protein (1-240AA)
Target	ATAD3B
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q5T9A4
Dilution Range	WB:1:500-5000, IHC-P:1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat(lane 1),...



Immunohistochemical staining of human ki...