


Product Datasheet

ADRBK1 antibody (orb352719)

Description	ADRBK1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Beta-adrenergic receptor kinase 1 protein (440-689AA)
Target	GRK2
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P25098
Dilution Range	WB:1:200-1000
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse brain tis...