


Product Datasheet

HSPB8 antibody (orb352728)

Description	HSPB8 antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Heat shock protein beta-8 protein (1-196AA)
Target	HSPB8
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9UJY1
Dilution Range	WB:1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of rat heart tissue...