



Product Datasheet

POTEH antibody (orb35279)

Description	Rabbit polyclonal antibody to POTEH. The protein is also known as A26C3 or ACTBL1 or POTE22, belongs...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Synthetic Peptide
Target	POTEH
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	60965
Uniprot ID	Q6S545
NCBI	NP_001129685.1
Dilution Range	WB: 1:1000, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Western blot analysis of in HeLa cell li...