



Product Datasheet

KLF3 antibody (orb352846)


Description	KLF3 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Recombinant Human Krueppel-like factor 3 protein (1-240AA)
Target	KLF3
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P57682
Dilution Range	IHC-P:1:20-1:100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human adrenal tissue.



Immunohistochemical staining of human kidney tissue.



Immunohistochemical staining of human brain tissue.