

Product Datasheet

N4BP2L2 antibody (FITC) (orb35285)

Description	N4BP2L2 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human NEDD4-binding protein 2-like 2 protein (1-300AA)
Target	N4BP2L2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q92802
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.