


Product Datasheet

OXCT1 antibody (orb352878)

Description	OXCT1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Succinyl-CoA:3-ketoacid coenzyme A transferase 1, mitochondrial protein (261-520AA)
Target	OXCT1
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P55809
Dilution Range	WB:1:500-5000
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse heart(Lan...