

Product Datasheet

C10orf46 antibody (orb35289)

Description	Rabbit polyclonal antibody to C10orf46. The protein is a member of the cullin family and is highly e...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	CACUL1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	41064
Uniprot ID	Q86Y37
NCBI	NP_722517.3
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

