



## Product Datasheet

### STK32A antibody (orb352906)


|                            |                                                                         |
|----------------------------|-------------------------------------------------------------------------|
| <b>Description</b>         | STK32A antibody                                                         |
| <b>Species/Host</b>        | Rabbit                                                                  |
| <b>Reactivity</b>          | Human                                                                   |
| <b>Conjugation</b>         | Unconjugated                                                            |
| <b>Tested Applications</b> | ELISA, IHC, WB                                                          |
| <b>Immunogen</b>           | Recombinant Human Serine/threonine-protein kinase 32A protein (1-166AA) |
| <b>Target</b>              | STK32A                                                                  |
| <b>Preservatives</b>       | PBS with 0.02% sodium azide, 50% glycerol, pH7.3.                       |
| <b>Form/Appearance</b>     | Liquid                                                                  |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.           |
| <b>Note</b>                | For research use only                                                   |
| <b>Isotype</b>             | IgG                                                                     |
| <b>Clonality</b>           | Polyclonal                                                              |
| <b>Uniprot ID</b>          | <b>Q8WU08</b>                                                           |
| <b>Dilution Range</b>      | WB:1:200-1000, IHC-P:1:20-1:200                                         |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                         |



Immunohistochemical staining of human br...



Immunohistochemical staining of human ad...



Western blot analysis of MCF-7 whole cel...