

Product Datasheet

NAA20 antibody (orb352908)

Description	NAA20 antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IP, WB
Immunogen	Recombinant Human N-alpha-acetyltransferase 20 protein (49-178AA)
Target	NAA20
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P61599
Dilution Range	WB:1:500-2000
Expiration Date	12 months from date of receipt.