



Product Datasheet

STAT5B antibody (orb352947)


Description	STAT5B antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ChIP, ELISA, IHC, WB
Immunogen	Recombinant Human Signal transducer and activator of transcription 5B protein (40-280AA)
Target	STAT5B
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P51692
Dilution Range	WB:1:500-2000, IHC-P:1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2(Lane 1), ...



Immunohistochemical staining of human te...



Immunoprecipitation analysis of MCF-7 ce...