

Product Datasheet

GABBR2 antibody (orb352960)

Description	GABBR2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human Gamma-aminobutyric acid type B receptor subunit 2 protein (60-310AA)
Target	GABBR2
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	075899
Expiration Date	12 months from date of receipt.