



## Product Datasheet

### SIRT3 antibody (orb352966)

|                            |                                                                                                 |
|----------------------------|-------------------------------------------------------------------------------------------------|
| <b>Description</b>         | SIRT3 antibody                                                                                  |
| <b>Species/Host</b>        | Rabbit                                                                                          |
| <b>Reactivity</b>          | Human                                                                                           |
| <b>Conjugation</b>         | Unconjugated                                                                                    |
| <b>Tested Applications</b> | ELISA, IHC, WB                                                                                  |
| <b>Immunogen</b>           | Recombinant Human NAD-dependent protein deacetylase sirtuin-3, mitochondrial protein (20-130AA) |
| <b>Target</b>              | SIRT3                                                                                           |
| <b>Preservatives</b>       | PBS with 0.02% sodium azide, 50% glycerol, pH7.3.                                               |
| <b>Form/Appearance</b>     | Liquid                                                                                          |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                   |
| <b>Note</b>                | For research use only                                                                           |
| <b>Isotype</b>             | IgG                                                                                             |
| <b>Clonality</b>           | Polyclonal                                                                                      |
| <b>Uniprot ID</b>          | <b>Q9NTG7</b>                                                                                   |
| <b>Dilution Range</b>      | WB:1:500-2000, IHC-P:1:20-1:200                                                                 |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                 |



Western blot analysis of NIH3T3(Lane 1),...



Immunohistochemical staining of breast c...